



## **FLOWRx Hub Manual**

### Accessible, affordable and scalable hub and spoke solutions

Our Manual FLOWRX Hub solution is designed for smaller groups of pharmacies looking to automate their original pack dispensing.

It allows pharmacies to free up enough capacity to start growing their businesses and spend more time on delivering private services and patient care. All solutions are scalable and additional stations and further automation can be added as your business grows.





In a FLOWRx Hub Manual, packs are scanned from the FLOW Tote to print a patient label, the label is affixed to the pack and scanned again. On this scan the RFID LED Smart Shelves guide the operator to sort the right medication into patientspecific baskets. You can then choose to automate the bagging or packing of medication by patient or leave this as a manual process too. Our software will then use barcode technology to sort patient medication bags into store specific totes.



# FLOW<sup>™</sup>H u b

#### At a glance.....

- Hub processing rate of up to **360 packs** per hour for a single line with one ACT/Pharmacist and two FTE
- Hub processing rate of **720 packs** per hour for a double line with one ACT/Pharmacist and five FTE.
- **Single Manual FLOWRx Hub:** Servicing somewhere between three and eight stores and processing 8,000 items per month with an 80/20% split for repeat prescriptions and walk-ins.
- **Double Manual FLOWRx Hub:** Servicing somewhere between eight and 17 stores and processing 8,000 items per month with an 80/20% split for repeat prescriptions and walk-ins.

### Layout and configuration

Manual Hub Single





